

Microdermabrasion+Electroporation

# AQUA TOUCH

The best combination for  
needle-free meso therapy

AquaTouch is the most dynamic skin treatment available today. It is in the only procedure featuring, the non-invasive exfoliation and delivery of skin-specific solutions, leaving patients with fresher-feeling and better-looking skin. AquaTouch is safe and painless, allowing patients to achieve optimum results on an accelerated basis without the complications and discomfort usually associated with invasive procedures. In addition, AquaTouch is unique because it does not use crystals or other abrasive particles that can irritate the skin.

## Benefits of the product

- Skin Refinement
- Skin Rejuvenation
- Reduction of Hyperpigmentation
- Fine Lines, Wrinkles & Sun Damage
- Improvement of Acne-prone skin & Puffy Eyes
- Demodex Removal
- Reduces wrinkles around the mouth



*No pain  
No side effect  
All skin type*

After the very first treatment, you will notice significant, visible results. Improve acne, hydrate dry, flaky skin and wrinkles and revitalize sun-damaged skin.

## How it works

The AquaTouch uses a closed loop vacuum system to exfoliate and apply topicals onto the skin. When its handpiece tip is brought in contact with the treatment area, the vacuum source gently pulls the tissue being treated through the hole in the tip and into contact with the abrasive, diamond-encrusted treatment head. As the handpiece is drawn across the skin, this abrasive surface exfoliates the tissue, applies a topical dermaceutical, and vacuums away the residuals into a sealed waste jar.

## Virtual mesotherapy without needles

Electrophoresis transfer the actives ingredients applied to the deeper layers using galvanic current. This allows performing mesotherapy without needles, enabling the realization of multiple treatments.

## Why is everybody raving about AquaTouch?

AquaTouch is the most dynamic skin treatment available today. It is the only procedure that provides both non-invasive exfoliation and delivery of skin-specific topical solutions under pressure, leaving patients with fresher-feeling and better-looking skin. AquaTouch is safe and painless without the complications and discomfort usually associated with crystal microdermabrasion, electroporation and other invasive treatments. In addition, AquaTouch is unique because it does not use crystals or other abrasive particles that can irritate the skin.

## Technical Specification

Power Input	AC220Vac 60 Hz
Maximum Power Consumption	100W
Classification (EN60601-1)	Class 1, Type BF
Mains Fuses	2 x T3A (5 x 20mm)
Size (height x width x depth)	400 x 310 x 900 mm
Weight	24kg
Max. Output Voltage	0~12V
Max. Output Current	0~300 $\mu$ A



## Accessories



**DAEYANG MEDICAL CO., LTD.**

Donghwa Medical Instrument Complex, 1656-1, Donghwa-ri, Munmak-eup, Wonju-si, Gangwon-do, 220-801 Korea  
Tel +82 33 747 4465~6 Fax +82 33 747 1133 E-mail daeyang@daeyangmed.com