

Non-Needle Fractional RF System

Liftron TherMatrix

Complete Skin Repair Renewal

No Surgery No Injection No Downtime

Features

- The Liftron features proprietary TherMatrix combined-energy technology of mono-polar and bi-polar radio frequency (RF).
- It's synergistically combined mono bipolar RF, improves safety for all skin types.
- Dual heating via mono & bi-polar RF energies raises target temperature while maintaining epidermal protection.
- The Liftron system consists of a single platform with a variety of treatment tips.
- Treat various body areas-such as the face, eyelids, lips, abdomen, arms, hands, thighs and buttocks.
- With multiple benefits-like skin tightening, contouring, body shaping and temporary improvement in the appearance of cellulite.
- All with one system-just change the tip.

Innovative Combined Treatment

- Non-invasive wrinkle treatment.
- Safe and effective treatments for the darker skin.
- Fast treatment, immediate visible results.
- A new category of fractional skin treatments.
- High dermal impact, low epidermal disruption for all skin types.
- Minimal downtime, proven efficacy, satisfied patients.



Feature of Non-Needle RF Tips

Dramatically improve texture, tone, laxity, wrinkles and skin irregularities with NFRF tips

1. Intelligently Designed tips to maximize treatment options and results
2. No Consumable or disposals
3. Multiple tip configurations

• Mono Fractional (MF) Tips

• Bi Fractional (BF) Tips



25 pins



100 pins



65 pins



100 pins



196 pins

Specification

- Model Name : Liftron TherMatrix
- Dimensions : 36(L) 32(D) 27(H) cm
- Weight : 5kg
- Model Name : Liftron DuoSupreme
- Dimensions : 47(W) 46(D) 100(H) cm
- Weight : 20kg
- Input : AC100/240V, 50/60Hz
- Max. Power Consumption : 20W
- Carrier Frequency : 1MHz
- Conduction Time : 50ms-300ms
- Run Time : 10-30m
- Applicator : Handprobe / Tip

